

## Datasheet: VPA00418

<b>Description:</b>	GOAT ANTI NOLC1
<b>Specificity:</b>	NOLC1
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Mouse  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid

#### Preparation

Goat polyclonal antibody purified by affinity chromatography

#### Buffer Solution

TRIS buffered saline

#### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)  
 1% Bovine Serum Albumin

#### Immunogen

Peptide sequence C-QVNSIKFDSE

---

**External Database****Links****UniProt:**[Q14978](#)   [Related reagents](#)**Entrez Gene:**[9221](#)   NOLC1   [Related reagents](#)

---

**Synonyms**

KIAA0035, NS5ATP13

---

**Specificity**

**Goat anti Human NOLC1 antibody** recognizes NOLC1, also known as nucleolar and coiled-body phosphoprotein 1, HCV NS5A-transactivated protein 13 and nucleolar protein p130.

Goat anti Human NOLC1 antibody detects a band of 130 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Anti NOLC1 detects a band of approximately 130 kDa in Jurkat cell lysate

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information**Material Safety Datasheet documentation #10058 available at:  
Antibody (10058): <https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf>

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)**North & South America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370758:200529'

Printed on 29 Aug 2021