

Datasheet: VPA00427

| | |
|----------------------|------------------------|
| Description: | RABBIT ANTI SRP14 |
| Specificity: | SRP14 |
| Format: | Purified |
| Product Type: | PrecisionAb Polyclonal |
| Isotype: | Polyclonal IgG |
| Quantity: | 100 µl |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ▪ | | | 1/1000 |

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

| | |
|---------------------------------|--|
| Target Species | Human |
| Product Form | Purified IgG - liquid. |
| Preparation | Rabbit polyclonal antibody purified by affinity chromatography. |
| Buffer Solution | Phosphate buffered saline. |
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃). |
| Immunogen | KLH conjugated synthetic peptide corresponding to amino acid 28-55 of human SRP14 |
| External Database Links | <p>UniProt: P37108 Related reagents</p> <p>Entrez Gene: 6727 SRP14 Related reagents</p> |

| | |
|--------------------------------------|---|
| Specificity | Rabbit anti Human SRP14 antibody recognizes SRP14 also known as signal recognition particle 14 kDa protein or 18 kDa Alu RNA-binding protein. Rabbit anti Human SRP14 antibody detects a band of 18 kDa. The antibody has been extensively validated for western blotting using whole cell lysates. |
| Western Blotting | Anti SRP14 detects a band of approximately 18 kDa in Ramos cell lysates. |
| Storage | Store undiluted at -20°C, avoiding repeated freeze thaw cycles. |
| Guarantee | 12 months from date of despatch. |
| Acknowledgements | PrecisionAb is a trademark of Bio-Rad Laboratories. |
| Health And Safety Information | Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf |
| Regulatory | For research purposes only. |

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M401639:220718'

Printed on 12 Sep 2022

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)