

## Datasheet: VPA00437

<b>Description:</b>	GOAT ANTI ARS2
<b>Specificity:</b>	ARS2
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Immunoprecipitation	▪			
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid.
<b>Preparation</b>	Goat polyclonal antibody purified by affinity chromatography.
<b>Buffer Solution</b>	Tris buffered saline.
<b>Preservative</b>	0.09% Sodium Azide (NaN <sub>3</sub> )
<b>Stabilisers</b>	0.5% Bovine Serum Albumin.
<b>Approx. Protein Concentrations</b>	IgG concentration 0.5 mg/ml
<b>Immunogen</b>	Peptide sequence C-DPRAIVEYRDLDPD
<b>External Database Links</b>	<b>UniProt:</b>

[Q9BXP5](#) [Related reagents](#)

**Entrez Gene:**

[51593](#) SRRT [Related reagents](#)

---

<b>Synonyms</b>	ARS2, ASR2
-----------------	------------

---

<b>Specificity</b>	<b>Goat anti Human ARS2 antibody</b> recognizes ARS2 also known as arsenate resistance protein 2, serrate RNA effector molecule homolog or SRRT.  Goat anti Human ARS2 antibody detects a band of 117 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.
--------------------	---

---

<b>Western Blotting</b>	Anti ARS2 detects a band of approximately 117 kDa in Jurkat cell lysates.
-------------------------	---

---

<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	---

---

<b>Guarantee</b>	12 months from date of despatch.
------------------	----------------------------------

---

<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories
-------------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10058 available at: Antibody (10058): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf</a>
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only.
-------------------	-----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M404435:220825'

Printed on 12 Sep 2022

---

© 2022 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)