

Datasheet: VPA00437

| | |
|----------------------|------------------------|
| Description: | GOAT ANTI ARS2 |
| Specificity: | ARS2 |
| Format: | Purified |
| Product Type: | PrecisionAb Polyclonal |
| Isotype: | Polyclonal IgG |
| Quantity: | 100 µl |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ▪ | | | 1/1000 |

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

| | |
|---------------------------------|---|
| Target Species | Human |
| Product Form | Purified IgG - liquid. |
| Preparation | Goat polyclonal antibody purified by affinity chromatography. |
| Buffer Solution | Tris buffered saline. |
| Preservative Stabilisers | 0.09% Sodium Azide (NaN ₃) 0.5% Bovine Serum Albumin. |
| Immunogen | Peptide sequence C-DPRAIVEYRDLDPD |
| External Database Links | UniProt: Q9BXP5 Related reagents Entrez Gene: |

| | |
|-----------------|------------|
| Synonyms | ARS2, ASR2 |
|-----------------|------------|

| | |
|--------------------|---|
| Specificity | Goat anti Human ARS2 antibody recognizes ARS2 also known as arsenate resistance protein 2, serrate RNA effector molecule homolog or SRRT. Goat anti Human ARS2 antibody detects a band of 117 kDa. The antibody has been extensively validated for western blotting using whole cell lysates. |
|--------------------|---|

| | |
|-------------------------|---|
| Western Blotting | Anti ARS2 detects a band of approximately 117 kDa in Jurkat cell lysates. |
|-------------------------|---|

| | |
|----------------|---|
| Storage | Store undiluted at -20°C, avoiding repeated freeze thaw cycles. |
|----------------|---|

| | |
|------------------|---------------------------------|
| Guarantee | 12 months from date of despatch |
|------------------|---------------------------------|

| | |
|-------------------------|--|
| Acknowledgements | PrecisionAb is a trademark of Bio-Rad Laboratories |
|-------------------------|--|

| | |
|--------------------------------------|--|
| Health And Safety Information | Material Safety Datasheet documentation #10058 available at: Antibody (10058): https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf |
|--------------------------------------|--|

| | |
|-------------------|-----------------------------|
| Regulatory | For research purposes only. |
|-------------------|-----------------------------|

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products.
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](https://www.bio-rad-antibodies.com)

'M370774:200529'

Printed on 12 Feb 2021