

Datasheet: VPA00459

BATCH NUMBER 172183

Description:	RABBIT ANTI SNRNP70
Specificity:	SNRNP70
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Learn about [how we validate our PrecisionAb range](#). Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid.

Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)
2% Sucrose.

Approx. Protein Concentrations	IgG concentration 0.5 mg/ml.
Immunogen	Synthetic peptide directed towards the N terminal region of human SNRP70
External Database Links	<p>UniProt: P08621 Related reagents</p> <p>Entrez Gene: 6625 SNRP70 Related reagents</p>
Synonyms	RNPU1Z, RPU1, SNRP70, U1AP1
Specificity	<p>Rabbit anti Human SNRP70 antibody recognizes SNRP70 also known as U1 small nuclear ribonucleoprotein 70 kDa, U1 snRNP 70 kDa.</p> <p>Rabbit anti Human SNRP70 antibody detects a band of 58 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti SNRP70 detects a band of approximately 58 kDa in Jurkat cell lysates.
Storage	This product is shipped at ambient temperature. Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Guarantee	12 months from date of despatch.
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10045 available at: https://www.bio-rad-antibodies.com/SDS/VPA00459
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M444277:250709'

Printed on 29 Jan 2026