

Datasheet: VPA00481

Description:	RABBIT ANTI RAP1A/B
Specificity:	RAP1A/B
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid.

Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)
2% Sucrose.

Immunogen

Synthetic peptide directed towards the N terminal region of human RAP1B

**External Database
Links**

UniProt:

[P61224](#) [Related reagents](#)

Entrez Gene:

[5908](#) RAP1B [Related reagents](#)

Specificity

Rabbit anti Human RAP1A / RAP1B antibody recognizes the Ras-related protein Rap-1a and Ras-related protein Rap-1b, also known as GTP-binding protein smg p21B, RAS-related protein RAP1B, Ras family small GTP binding protein RAP1B, small GTP binding protein.

RAP1A/RAP1B encodes members of the RAS-like small GTP-binding protein superfamily. Members of this family regulate multiple cellular processes including cell adhesion and growth and differentiation. The RAP1A protein regulates signaling pathways that affect cell proliferation and adhesion, and may play a role in tumor malignancy. The RAP1B protein localizes to cellular membranes and has been shown to regulate integrin-mediated cell signaling. RAP1B also plays a role in regulating outside-in signaling in platelets. Alternate splicing of RAP1A and RAP1B have been known to produce multiple transcript variants. Pseudogenes of RAP1A have been defined on chromosomes 14 and 17 (provided by RefSeq, May 2014). Pseudogenes of RAP1B are found on chromosomes 3, 5, 6 and 9. (provided by RefSeq, Oct 2011).

Rabbit anti Human RAP1A/ RAP1B antibody detects a band of 26 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti RAP1A/B detects a band of approximately 26 kDa in HEK293 cell lysates.

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch.

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety
Information**

Material Safety Datasheet documentation #10045 available at:
<https://www.bio-rad-antibodies.com/SDS/VPA00481>
Antibody (10045)

Regulatory

For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M402303:220718'

Printed on 13 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)