

## Datasheet: VPA00503

<b>Description:</b>	RABBIT ANTI RPL12
<b>Specificity:</b>	RPL12
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Mouse  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid.

#### Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

#### Buffer Solution

Phosphate buffered saline.

#### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)  
 2% Sucrose.

#### Immunogen

Synthetic peptide directed towards the N terminal region of human RPL12.

---

**External Database****Links****UniProt:**

[P30050](#)   [Related reagents](#)

**Entrez Gene:**

[6136](#)   RPL12   [Related reagents](#)

---

**Specificity**

**Rabbit anti Human RPL12 antibody** recognizes RPL12 also known as 60S ribosomal protein L12.

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. RPL12 gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L11P family of ribosomal proteins. It is located in the cytoplasm. The protein binds directly to the 26S rRNA. RPL12 is co-transcribed with the U65 snoRNA, which is located in its fourth intron. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of RPL12 dispersed through the genome (provided by RefSeq, Jul 2008).

Rabbit anti Human RPL12 antibody detects a band of 21 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Anti RPL12 detects a band of approximately 21 kDa in HEK293 cell lysates.

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information**

Material Safety Datasheet documentation #10045 available at:  
<https://www.bio-rad-antibodies.com/SDS/VPA00503>  
Antibody (10045)

---

**Regulatory**

For research purposes only.

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South**   Tel: +1 800 265 7376

**America**   Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370823:200529'

Printed on 13 Aug 2023

---

