

## Datasheet: VPA00514

<b>Description:</b>	GOAT ANTI HRS / HGS
<b>Specificity:</b>	HRS / HGS
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Mouse

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid.

#### Preparation

Goat polyclonal antibody purified by affinity chromatography.

#### Buffer Solution

TRIS buffered saline.

#### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)  
0.5% Bovine Serum Albumin.

#### Immunogen

Peptide sequence C-NRNYWEKKQEEARK

**External Database  
Links**

**UniProt:**

[O14964](#)   [Related reagents](#)

**Entrez Gene:**

[9146](#)   HGS   [Related reagents](#)

---

**Synonyms**                    HRS

---

**Specificity**                    **Goat anti Human HRS / HGS antibody** recognizes HRS / HGS, also known as human growth factor-regulated tyrosine kinase substrate or protein pp110.

The protein encoded by HGS regulates endosomal sorting and plays a critical role in the recycling and degradation of membrane receptors. The encoded protein sorts monoubiquitinated membrane proteins into the multivesicular body, targeting these proteins for lysosome-dependent degradation (provided by RefSeq, Dec 2010).

Goat anti Human HRS / HGS antibody detects a band of 110 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**            Anti HRS / HGS detects a band of approximately 110 kDa in HEK293 cell lysates.

---

**Storage**                        Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee**                    12 months from date of despatch.

---

**Acknowledgements**        PrecisionAb is a trademark of Bio-Rad Laboratories.

---

**Health And Safety  
Information**                    Material Safety Datasheet documentation #10058 available at:  
Antibody (10058): <https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf>

---

**Regulatory**                    For research purposes only.

---

## Related Products

### Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

**North & South**    Tel: +1 800 265 7376

**America**            Fax: +1 919 878 3751

                          Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

                          Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

                          Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M402430:220719'

**Printed on 19 Jul 2022**