

Datasheet: VPA00529

| | |
|----------------------|------------------------|
| Description: | RABBIT ANTI ACTN3 |
| Specificity: | ACTN3 |
| Format: | Purified |
| Product Type: | PrecisionAb Polyclonal |
| Isotype: | Polyclonal IgG |
| Quantity: | 100 µl |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ▪ | | | 1/1000 |

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse, Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid.

Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

Buffer Solution

Phosphate buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)
2% Sucrose.

Immunogen

Synthetic peptide directed towards the N terminal region of human ACTN3

**External Database
Links**

UniProt:

[Q08043](#) [Related reagents](#)

Entrez Gene:

[89](#) ACTN3 [Related reagents](#)

Specificity

Rabbit anti Human ACTN3 antibody recognizes ACTN3 also known as alpha-actinin-3, F-actin cross-linking protein, alpha-actinin skeletal muscle isoform 3.

ACTN3 gene encodes a member of the alpha-actin binding protein gene family. The encoded protein is primarily expressed in skeletal muscle and functions as a structural component of sarcomeric Z line. This protein is involved in crosslinking actin containing thin filaments. An allelic polymorphism in ACTN3 results in both coding and non-coding variants. The reference genome contains the non-coding allele. The non-functional allele of ACTN3 is associated with elite athlete status (provided by RefSeq, Apr 2012).

Rabbit anti Human ACTN3 antibody detects a band of 103 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti ACTN3 detects a band of approximately 103 kDa in NIH/3T3 cell lysates.

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch.

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety
Information**

Material Safety Datasheet documentation #10045 available at:
<https://www.bio-rad-antibodies.com/SDS/VPA00529>
Antibody (10045)

Regulatory

For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M402492:220719'

Printed on 13 Aug 2023