

## Datasheet: VPA00532

<b>Description:</b>	RABBIT ANTI BETA-CENTRACTIN
<b>Specificity:</b>	BETA-CENTRACTIN
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Mouse <b>N.B.</b> Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.
<b>Product Form</b>	Purified IgG - liquid.
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline.
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ) 2% Sucrose.
<b>Immunogen</b>	Synthetic peptide directed towards the C terminal region of human beta-centractin

**External Database  
Links**

**UniProt:**

[P42025](#)

[Related reagents](#)

**Entrez Gene:**

[10120](#)

ACTR1B

[Related reagents](#)

**Synonyms**

CTRN2

**Specificity**

**Rabbit anti Human beta-centractin antibody** recognizes beta-centractin also known as actin-related protein 1B or ACTR1B.

The ACTR1B gene encodes a 42.3 kDa subunit of dynactin, a macromolecular complex consisting of 10 subunits ranging in size from 22 to 150 kDa. Dynactin binds to both microtubules and cytoplasmic dynein and is involved in a diverse array of cellular functions, including ER-to-Golgi transport, the centripetal movement of lysosomes and endosomes, spindle formation, chromosome movement, nuclear positioning, and axonogenesis. This subunit, like ACTR1A, is an actin-related protein. These two proteins, which are of equal length and share 90% amino acid identity, are present in a constant ratio of approximately 1:15 in the dynactin complex (provided by RefSeq, Aug 2008).

Rabbit anti Human beta-centractin antibody detects a band of 42 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

**Western Blotting**

Anti beta-centractin detects a band of approximately 42 kDa in HeLa cell lysates.

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

**Guarantee**

12 months from date of despatch.

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety  
Information**

Material Safety Datasheet documentation #10045 available at:  
<https://www.bio-rad-antibodies.com/SDS/VPA00532>  
Antibody (10045)

**Regulatory**

For research purposes only.

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South  
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M402504:220719'

