

## Datasheet: VPA00542

<b>Description:</b>	RABBIT ANTI SEL1L
<b>Specificity:</b>	SEL1L
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Mouse  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid.

#### Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

#### Buffer Solution

Phosphate buffered saline.

#### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>).

#### Immunogen

KLH conjugated synthetic peptide corresponding to amino acid 156-182 of human SEL1L

<b>External Database Links</b>	<b>UniProt:</b> <a href="#">Q9UBV2</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">6400</a> SEL1L <a href="#">Related reagents</a>
<b>Synonyms</b>	TSA305
<b>Specificity</b>	<p><b>Rabbit anti Human SEL1L antibody</b> recognizes SEL1L, also known as suppressor of lin 12 or suppressor of lin-12-like protein 1.</p> <p>The protein encoded by SEL1L gene is part of a protein complex required for the retrotranslocation or dislocation of misfolded proteins from the endoplasmic reticulum lumen to the cytosol, where they are degraded by the proteasome in a ubiquitin-dependent manner. Alternatively spliced transcript variants encoding different isoforms have been found for SEL1L (provided by RefSeq, Oct 2011).</p> <p>Rabbit anti Human SEL1L antibody detects a band of 89 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti SEL1L detects a band of approximately 89 kDa in HeLa cell lysates.
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 available at: Antibody (10040): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a>
<b>Regulatory</b>	For research purposes only.

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M370850:200529'

Printed on 07 May 2021