

Datasheet: VPA00548

BATCH NUMBER 160712

Description:	GOAT ANTI CYCLOPHILIN A
Specificity:	CYCLOPHILIN A
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid.

Preparation

Goat polyclonal antibody purified by affinity chromatography.

Buffer Solution

TRIS buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)

0.5% Bovine Serum Albumin.

Immunogen Peptide sequence C-ERFGSRNGKTSKK

External Database

Links

UniProt:

[P62937](#) [Related reagents](#)

Entrez Gene:

[5478](#) PPIA [Related reagents](#)

Synonyms

CYPA

Specificity

Goat anti Human cyclophilin A antibody recognizes cyclophilin A also known as PPIase A, PPIA, T cell cyclophilin, cyclosporin A-binding protein, epididymis secretory sperm binding protein Li 69p, peptidyl-prolyl cis-trans isomerase A or rotamase A.

PPIA gene encodes a member of the peptidyl-prolyl cis-trans isomerase (PPIase) family. PPIases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. The encoded protein is a cyclosporin binding-protein and may play a role in cyclosporin A-mediated immunosuppression. Cyclophilin A can also interact with several HIV proteins, including p55 gag, Vpr, and capsid protein, and has been shown to be necessary for the formation of infectious HIV virions ([Braaten et al. 1996](#)). Multiple pseudogenes that map to different chromosomes have been reported (provided by RefSeq, Jul 2008).

Goat anti Human cyclophilin A antibody detects a band of 18 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti cyclophilin A detects a band of approximately 18 kDa in A549 cell lysates.

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories

Health And Safety Information

Material Safety Datasheet documentation #10058 available at: <https://www.bio-rad-antibodies.com/SDS/VPA00548>
Antibody (10058)

Regulatory

For research purposes only.

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M370856:200529'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)