

Datasheet: VPA00552

| | |
|----------------------|--|
| Description: | GOAT ANTI RE1-SILENCING TRANSCRIPTION FACTOR |
| Specificity: | RE1-SILENCING TRANSCRIPTION FACTOR |
| Format: | Purified |
| Product Type: | PrecisionAb Polyclonal |
| Isotype: | Polyclonal IgG |
| Quantity: | 100 µl |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ▪ | | | 1/1000 |

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid.

Preparation

Goat polyclonal antibody purified by affinity chromatography.

Buffer Solution

TRIS buffered saline.

Preservative Stabilisers

0.09% Sodium Azide (NaN₃)
0.5% Bovine Serum Albumin.

Immunogen

Peptide sequence C-PETPGAEDKKGSSK

**External Database
Links**

UniProt:

[Q13127](#) [Related reagents](#)

Entrez Gene:

[5978](#) REST [Related reagents](#)

Synonyms

NRSF, XBR

Specificity

Goat anti Human RE1-silencing transcription factor antibody recognizes the RE1-silencing transcription factor also known as REST, neural-restrictive silencer factor or X2 box repressor.

The REST gene encodes a transcriptional repressor that represses neuronal genes in non-neuronal tissues. It is a member of the Kruppel-type zinc finger transcription factor family. It represses transcription by binding a DNA sequence element called the neuron-restrictive silencer element. The protein is also found in undifferentiated neuronal progenitor cells and it is thought that this repressor may act as a master negative regulator of neurogenesis. Alternatively spliced transcript variants have been described (provided by RefSeq, Jul 2010).

Goat anti Human RE1-silencing transcription factor antibody detects a band of 121 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti RE1-silencing transcription factor detects a band of approximately 121 kDa in Raji cell lysates.

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee

12 months from date of despatch.

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories.

**Health And Safety
Information**

Material Safety Datasheet documentation #10058 available at:
<https://www.bio-rad-antibodies.com/SDS/VPA00552>
Antibody (10058)

Regulatory

For research purposes only.

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

**North & South
America** Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M402595:220719'

