

Datasheet: VPA00564

Description:	RABBIT ANTI SPARC
Specificity:	SPARC
Format:	Purified
Product Type:	PrecisionAb™ Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid.
Preparation	Rabbit polyclonal antibody purified by affinity chromatography.
Buffer Solution	Phosphate buffered saline.
Preservative	0.09% Sodium Azide (NaN ₃)
Stabilisers	2% Sucrose.
Immunogen	Synthetic peptide directed towards the C terminal region of human SPARC
External Database Links	<p>UniProt: P09486 Related reagents</p> <p>Entrez Gene: 6678 SPARC Related reagents</p>
Synonyms	ON

Specificity **Rabbit anti Human SPARC antibody** recognizes SPARC, also known as BM-40, basement-membrane protein 40, osteonectin or secreted protein acidic and rich in cysteine.

The SPARC gene encodes a cysteine-rich acidic matrix-associated protein. The encoded protein is

required for the collagen in bone to become calcified but is also involved in extracellular matrix synthesis and promotion of changes to cell shape. The gene product has been associated with tumor suppression but has also been correlated with metastasis based on changes to cell shape which can promote tumor cell invasion (provided by RefSeq, Dec 2011). Mutation of the SPARC gene causes presentation of Osteogenesis imperfecta 17 ([OI17](#)), an autosomal recessive form of osteogenesis imperfecta, a connective tissue disorder characterized by low bone mass, bone fragility and susceptibility to fractures after minimal trauma ([Mendoza-Londono et al. 2015](#)).

Rabbit anti Human SPARC antibody detects a band of 35 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Anti SPARC detects a band of approximately 35 kDa in HEK293 cell lysates.

Instructions For Use Please refer to the [PrecisionAb western blotting protocol](#).

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee As supplied, 12 months from date of despatch.

Acknowledgements PrecisionAb™ is a trademark of Bio-Rad Laboratories.

Health And Safety Information Material Safety Datasheet documentation #10045 available at: Antibody (10045): <https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf>

Regulatory For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M338542:181217'

Printed on 20 May 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)