

## Datasheet: VPA00570KT

<b>Description:</b>	PIAS3 ANTIBODY WITH CONTROL LYSATE
<b>Specificity:</b>	PIAS3
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	2 Westerns

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**PrecisionAb antibodies have been extensively [validated for the western blot application](#).** The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Mouse  
**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid.

#### Preparation

20µl Rabbit Ig fraction prepared by ammonium sulphate precipitation.

#### Buffer Solution

Phosphate buffered saline.

#### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>).

#### Immunogen

KLH conjugated synthetic peptide corresponding to amino acid 95-126 of human PIAS3

<b>External Database Links</b>	<b>UniProt:</b> <a href="#">Q9Y6X2</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">10401</a> PIAS3 <a href="#">Related reagents</a>
<b>Specificity</b>	<p><b>Rabbit anti Human PIAS3 antibody</b> recognizes PIAS3, also known as E3 SUMO-protein ligase PIAS3, protein inhibitor of activated STAT protein 3 or "zinc finger, MIZ-type containing 5".</p> <p>The PIAS3 gene encodes a member of the PIAS (protein inhibitor of activated STAT (signal transducer and activator of transcription) family of transcriptional modulators. The protein functions as a SUMO (small ubiquitin-like modifier)-E3 ligase which catalyzes the covalent attachment of a SUMO protein to specific target substrates. It directly binds to several transcription factors and either blocks or enhances their activity. Alternatively spliced transcript variants of PIAS3 have been identified, but the full-length nature of some of these variants has not been determined (provided by RefSeq, Jul 2008).</p> <p>Rabbit anti Human PIAS3 antibody detects a band of 75 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti PIAS3 detects a band of approximately 75 kDa in MCF-7 cell lysates.
<b>Instructions For Use</b>	Please refer to the PrecisionAb western blotting protocol.
<b>Lysate Composition</b>	400µg MCF7 lysate lyophilized in RIPA buffer.
<b>Lysate Reconstitution</b>	<ul style="list-style-type: none"> <li>- If using DDT reconstitute the lyophilized lysate with 190µl DI H<sub>2</sub>O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.</li> <li>- If using BME reconstitute the lyophilized lysate with 180µl DI H<sub>2</sub>O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.</li> </ul> <p>Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.</p>
<b>Storage</b>	<p>Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.</p> <p>Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.</p>
<b>Guarantee</b>	As supplied, 12 months from date of despatch.
<b>Acknowledgements</b>	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10040 #10561 available at: Antibody (10040): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</a> Lysate Material (10561): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf</a>

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

From March 15, 2021, we will no longer supply printed datasheets with our products.  
Look out for updates on how to access your digital version at [bio-rad-antibodies.com](http://bio-rad-antibodies.com)

'M351322:190318'

**Printed on 10 Feb 2021**

---