

Datasheet: VPA00578

**BATCH NUMBER 160811**

<b>Description:</b>	RABBIT ANTI JAB1
<b>Specificity:</b>	JAB1
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid.

### Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

### Buffer Solution

Phosphate buffered saline.

### Preservative Stabilisers

0.09% Sodium Azide (NaN<sub>3</sub>)  
2% Sucrose.

**Immunogen** Synthetic peptide directed towards the N terminal region of human JAB1

---

**External Database**

**Links**

**UniProt:**

[Q92905](#)    [Related reagents](#)

**Entrez Gene:**

[10987](#)    COPS5    [Related reagents](#)

---

**Synonyms**

CSN5, JAB1

---

**Specificity**

**Rabbit anti Human JAB1 antibody** recognizes JAB1 also known as COP9 signalosome complex subunit 5, COPS5, 38 kDa Mov34 homolog, COP9 constitutive photomorphogenic homolog subunit 5, jun activation domain-binding protein 1 or signalosome subunit 5.

The protein encoded by the COPS5 gene is one of the eight subunits of COP9 signalosome, a highly conserved protein complex that functions as an important regulator in multiple signaling pathways. The structure and function of COP9 signalosome is similar to that of the 19S regulatory particle of 26S proteasome. COP9 signalosome has been shown to interact with SCF-type E3 ubiquitin ligases and act as a positive regulator of E3 ubiquitin ligases. This protein is reported to be involved in the degradation of cyclin-dependent kinase inhibitor CDKN1B/p27Kip1. It is also known to be an coactivator that increases the specificity of JUP1 transcription factors (provided by RefSeq, Jul 2008).

Rabbit anti Human JAB1 antibody detects a band of 37 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Anti JAB1 detects a band of approximately 37 kDa in HeLa cell lysates.

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information**

Material Safety Datasheet documentation #10045 available at:  
Antibody (10045): <https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf>

---

**Regulatory**

For research purposes only.

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South**    Tel: +1 800 265 7376

**America**    Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M370883:200529'

**Printed on 08 Feb 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)