

Datasheet: VPA00591KT

Description:	PLSCR1 ANTIBODY WITH CONTROL LYSATE
Specificity:	PLSCR1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	2 Westerns

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid.
Preparation	20µl Rabbit polyclonal antibody purified by affinity chromatography.
Buffer Solution	Phosphate buffered saline.
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 2% Sucrose.
Immunogen	Synthetic peptide directed towards the N terminal region of human PLSCR1

External Database Links

UniProt:

[O15162](#)

[Related reagents](#)

Entrez Gene:

[5359](#)

PLSCR1

[Related reagents](#)

Specificity	<p>Rabbit anti Human PLSCR1 antibody recognizes PLSCR1 also known as phospholipid scramblase 1, PL scramblase 1, calcium dependent phospholipid scramblase 1, erythrocyte phospholipid scramblase or mTRA1b. Phospholipid scramblase 1 is a 318 amino acid ~35 kDa single pass type II transmembrane protein</p> <p>Rabbit anti Human PLSCR1 antibody detects a band of 35 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti PLSCR1 detects a band of approximately 35 kDa in A431 cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol .
Lysate Composition	400µg A431 lysate lyophilized in RIPA buffer.
Lysate Reconstitution	<p>- If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.</p> <p>- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.</p> <p>Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.</p>
Storage	<p>Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.</p> <p>Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.</p>
Guarantee	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	<p>Material Safety Datasheet documentation #10045 #10561 available at:</p> <p>Antibody (10045): https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf</p> <p>Lysate Material (10561): https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf</p>
Regulatory	For research purposes only.

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

From March 15, 2021, we will no longer supply printed datasheets with our products. Look out for updates on how to access your digital version at bio-rad-antibodies.com

'M351355:190318'

Printed on 10 Feb 2021

