

## Datasheet: VPA00613

<b>Description:</b>	RABBIT ANTI RXRB
<b>Specificity:</b>	RXRB
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide 2% Sucrose
<b>Immunogen</b>	Synthetic peptide encompassing part of the N terminal region of human RXRB
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P28702</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b>

---

<b>Synonyms</b>	NR2B2
-----------------	-------

---

<b>Specificity</b>	<p><b>Rabbit anti Human RXRB antibody</b> recognizes the retinoic acid receptor RXR-beta, also known as MHC class I promoter binding protein or nuclear receptor subfamily 2 group B member 2.</p> <p>The RXRB gene encodes a member of the retinoid X receptor (RXR) family of nuclear receptors which are involved in mediating the effects of retinoic acid (RA). The encoded protein forms homodimers with the retinoic acid, thyroid hormone, and vitamin D receptors, increasing both DNA binding and transcriptional function on their respective response elements. RXRB lies within the major histocompatibility complex (MHC) class II region on chromosome 6. Alternatively spliced transcript variants encoding multiple isoforms have been observed for RXRB (provided by RefSeq, Jul 2012).</p> <p>Rabbit anti Human RXRB antibody detects a band of 70 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
--------------------	---

---

<b>Western Blotting</b>	Rabbit anti RXRB detects a band of approximately 70 kDa in K562 cell lysates
-------------------------	--

---

<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	---

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories
-------------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10045 available at: Antibody (10045): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf</a>
--------------------------------------	--

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370913:200529'

Printed on 29 Aug 2021