

Datasheet: VPA00621

**BATCH NUMBER 161121**

<b>Description:</b>	RABBIT ANTI CHMP2A
<b>Specificity:</b>	CHMP2A
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Rabbit polyclonal antibody purified by affinity chromatography

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide  
2% Sucrose

**Immunogen** Synthetic peptide encompassing part of the middle region of human CHMP2A

---

**External Database**

**Links**

**UniProt:**

[O43633](#)

[Related reagents](#)

**Entrez Gene:**

[27243](#)

CHMP2A

[Related reagents](#)

---

**Synonyms**

BC2, CHMP2

---

**Specificity**

**Rabbit anti Human CHMP2A antibody** recognizes the charged multivesicular body protein 2a, also known as VPS2 homolog A, chromatin modifying protein 2A, putative breast adenocarcinoma marker BC-2 (32kD) or vacuolar protein sorting-associated protein(vps) 2-1.

CHMP2A belongs to the chromatin-modifying protein/charged multivesicular body protein (CHMP) family. These proteins are components of ESCRT-III (endosomal sorting complex required for transport III), a complex involved in degradation of surface receptor proteins and formation of endocytic multivesicular bodies (MVBs). Some CHMPs have both nuclear and cytoplasmic/vesicular distributions, and one such CHMP, CHMP1A (MIM 164010), is required for both MVB formation and regulation of cell cycle progression ([Tsang \*et al.\* 2006](#)).

Rabbit anti Human CHMP2A antibody detects a band of 25 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Rabbit anti CHMP2A detects a band of approximately 25 kDa in Raji cell lysates

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information**

Material Safety Datasheet documentation #10045 available at:

<https://www.bio-rad-antibodies.com/SDS/VPA00621>

Antibody (10045)

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M370921:200529'

**Printed on 13 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)