

## Datasheet: VPA00622

<b>Description:</b>	RABBIT ANTI CDC73
<b>Specificity:</b>	CDC73
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

#### Target Species

Human

#### Species Cross Reactivity

Reacts with: Mouse, Rat

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

#### Product Form

Purified IgG - liquid

#### Preparation

Rabbit polyclonal antibody purified by affinity chromatography.

#### Buffer Solution

Phosphate buffered saline.

#### Preservative Stabilisers

0.09% Sodium Azide  
2% Sucrose

#### Immunogen

Synthetic peptide encompassing part of the middle region of human CDC73

**External Database  
Links**

**UniProt:**

[Q6P1J9](#) [Related reagents](#)

**Entrez Gene:**

[79577](#) CDC73 [Related reagents](#)

---

**Synonyms**

C1orf28, HRPT2

---

**Specificity**

**Rabbit anti Human CDC73 antibody** recognizes parafibromin, also known as also known as Paf1/RNA polymerase II complex component cell division cycle 73, cell division cycle protein 73 homolog, HRPT2 or hyperparathyroidism 2 protein.

CDC73 gene encodes a tumor suppressor that is involved in transcriptional and post-transcriptional control pathways. The protein is a component of the the PAF protein complex, which associates with the RNA polymerase II subunit POLR2A and with a histone methyltransferase complex. This protein appears to facilitate the association of 3' mRNA processing factors with actively-transcribed chromatin. Mutations in CDC73 have been linked to hyperparathyroidism-jaw tumor syndrome, familial isolated hyperparathyroidism, and parathyroid carcinoma (provided by RefSeq, Jul 2009).

Rabbit anti Human CDC73 antibody detects a band of 60 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Rabbit anti CDC73 detects a band of approximately 60 kDa in NIH/3T3 cell lysates.

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee**

12 months from date of despatch.

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories.

---

**Health And Safety  
Information**

Material Safety Datasheet documentation #10045 available at:  
Antibody (10045): <https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf>

---

**Regulatory**

For research purposes only.

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South**

Tel: +1 800 265 7376

**America**

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M402907:220720'

Printed on 07 Feb 2023

