

Datasheet: VPA00639KT

| | |
|----------------------|----------------------------------|
| Description: | CAD ANTIBODY WITH CONTROL LYSATE |
| Specificity: | CAD |
| Format: | Purified |
| Product Type: | PrecisionAb™ Polyclonal |
| Isotype: | Polyclonal IgG |
| Quantity: | 2 Westerns |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ▪ | | | 1/1000 |

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

| | |
|--------------------------------|--|
| Target Species | Human |
| Product Form | Purified IgG - liquid |
| Preparation | 20µl Rabbit polyclonal antibody purified by affinity chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative | 0.09% Sodium Azide |
| Stabilisers | 2% Sucrose |
| Immunogen | Synthetic peptide encompassing part of the middle region of human CAD |
| External Database Links | <p>UniProt: O76075 Related reagents</p> <p>Entrez Gene: 1677 DFFB Related reagents</p> |
| Synonyms | CAD, DFF2, DFF40 |
| Specificity | Rabbit anti Human CAD antibody recognizes the DNA fragmentation factor subunit beta, also known as DNA fragmentation factor 40 kDa subunit, caspase-activated DNase, DFF40, CPAN or caspase-activated nuclease. |

Apoptosis is a cell death process that removes toxic and/or useless cells during mammalian development. The apoptotic process is accompanied by shrinkage and fragmentation of the cells and nuclei and degradation of the chromosomal DNA into nucleosomal units. DNA fragmentation factor (DFF) is a heterodimeric protein of 40-kD (DFFB) and 45-kD (DFFA) subunits. DFFA is the substrate for caspase-3 and triggers DNA fragmentation during apoptosis. DFF becomes activated when DFFA is cleaved by caspase-3. The cleaved fragments of DFFA dissociate from DFFB, the active component of DFF. DFFB has been found to trigger both DNA fragmentation and chromatin condensation during apoptosis. Alternatively spliced transcript variants encoding distinct isoforms have been found for DFFB but the biological validity of some of these variants has not been determined (provided by RefSeq, Sep 2013).

Rabbit anti Human CAD antibody detects a band of 46 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Rabbit anti CAD detects a band of approximately 46 kDa in HeLa cell lysates

Instructions For Use Please refer to the [PrecisionAb western blotting protocol](#). For additional information on secondary antibody dilution and exposure time see product web page.

Lysate Composition 400µg HeLa lysate lyophilized in RIPA buffer.

Lysate Reconstitution - If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.
- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.
Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.

Storage Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.

Shelf Life As supplied, 12 months from date of despatch

Acknowledgements PrecisionAb™ is a trademark of Bio-Rad Laboratories

Health And Safety Information Material Safety Datasheet documentation #10045 #10561 available at:
Antibody (10045): <https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf>
Lysate Material (10561): <https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

