

## Datasheet: VPA00645

<b>Description:</b>	RABBIT ANTI NESTIN
<b>Specificity:</b>	NESTIN
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide 2% Sucrose
<b>Immunogen</b>	Synthetic peptide encompassing part of the middle region of human nestin

#### External Database Links

**UniProt:**  
[P48681](#)    [Related reagents](#)

**Entrez Gene:**

---

**Specificity** Rabbit anti Human nestin antibody recognizes nestin.

The NES gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells (provided by RefSeq, Sep 2011).

Rabbit anti Human nestin antibody detects a band of 290 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting** Rabbit anti nestin detects a band of approximately 290 kDa in NTERA-2 cell lysates

---

**Storage** Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

---

**Guarantee** 12 months from date of despatch

---

**Acknowledgements** PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information** Material Safety Datasheet documentation #10045 available at: Antibody (10045): <https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf>

---

**Regulatory** For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M370941:200529'

Printed on 22 Mar 2021