

## Datasheet: VPA00657

<b>Description:</b>	RABBIT ANTI RPRD2
<b>Specificity:</b>	RPRD2
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography.
<b>Buffer Solution</b>	Phosphate buffered saline.
<b>Preservative</b>	0.09% Sodium Azide
<b>Stabilisers</b>	2% Sucrose
<b>Immunogen</b>	Synthetic peptide encompassing part of the middle region of human RPRD2

#### External Database Links

##### UniProt:

[Q5VT52](#)

[Related reagents](#)

##### Entrez Gene:

[23248](#)

RPRD2

[Related reagents](#)

<b>Synonyms</b>	KIAA0460
<b>Specificity</b>	<b>Rabbit anti Human RPRD2 antibody</b> recognizes the regulation of nuclear pre-mRNA domain-containing protein 2, also known as regulation of nuclear pre-mRNA domain-containing protein 2. This protein associates to RNA Pol II to promote snRNA transcription.  Rabbit anti Human RPRD2 antibody detects a band of 220 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.
<b>Western Blotting</b>	Rabbit anti RPRD2 detects a band of approximately 220 kDa in Jurkat cell lysates.
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
<b>Guarantee</b>	12 months from date of despatch.
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10045 available at: <a href="https://www.bio-rad-antibodies.com/SDS/VPA00657">https://www.bio-rad-antibodies.com/SDS/VPA00657</a> Antibody (10045)
<b>Regulatory</b>	For research purposes only.

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M416402:230301'

Printed on 13 Aug 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)