

Datasheet: VPA00665

**BATCH NUMBER 161118**

<b>Description:</b>	RABBIT ANTI APOLIPOPROTEIN E
<b>Specificity:</b>	APOLIPOPROTEIN E
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.01% Thiomersal 1% Bovine Serum Albumin <50% Glycerol
<b>Immunogen</b>	Carrier-protein conjugated synthetic peptide encompassing part of the C terminal region of human apolipoprotein E
<b>External Database Links</b>	<b>UniProt:</b>

**Entrez Gene:**[348](#)

APOE

[Related reagents](#)**Specificity**

**Rabbit anti Human apolipoprotein E antibody** recognizes apolipoprotein E, also known as apolipoprotein E3, APOE, APO-E, or LPG.

Chylomicron remnants and very low density lipoprotein (VLDL) remnants are rapidly removed from the circulation by receptor-mediated endocytosis in the liver. Apolipoprotein E, a main apoprotein of the chylomicron, binds to a specific receptor on liver cells and peripheral cells. ApoE is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. The APOE gene is mapped to chromosome 19 in a cluster with APOC1 and APOC2. Defects in apolipoprotein E result in familial dysbetalipoproteinemia, or type III hyperlipoproteinemia (HLP III), in which increased plasma cholesterol and triglycerides are the consequence of impaired clearance of chylomicron and VLDL remnants (provided by RefSeq, Jul 2008).

Rabbit anti Human apolipoprotein E antibody detects a band of 35 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

**Western Blotting**

Rabbit anti apolipoprotein E detects a band of approximately 35 kDa in HepG2 cell lysates

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

**Guarantee**

12 months from date of despatch

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories

**Health And Safety Information**

Material Safety Datasheet documentation #10096 available at:  
<https://www.bio-rad-antibodies.com/SDS/VPA00665>  
Antibody (10096)

**Regulatory**

For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370956:200529'

Printed on 19 Jan 2024