

Datasheet: VPA00707

Description:	RABBIT ANTI PEROXIREDOXIN 6
Specificity:	PEROXIREDOXIN 6
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by affinity chromatography on Protein A
Buffer Solution	HEPES buffered saline
Preservative Stabilisers	0.09% Sodium Azide 1% Bovine Serum Albumin <50% Glycerol
Immunogen	<i>E. coli</i> -derived recombinant human peroxiredoxin 6
External Database Links	UniProt: P30041 Related reagents

Entrez Gene:

[9588](#) PRDX6 [Related reagents](#)

Synonyms AOP2, KIAA0106

Specificity **Rabbit anti Human peroxiredoxin 6 antibody** recognizes peroxiredoxin-6, also known as 1-Cys PRX, 1-Cys peroxiredoxin, acidic calcium-independent phospholipase A2, PRDX6, or non-selenium glutathione peroxidase.

The protein encoded by PRDX6 gene is a member of the thiol-specific antioxidant protein family. This protein is a bifunctional enzyme with two distinct active sites. It is involved in redox regulation of the cell; it can reduce H₂O₂ and short chain organic, fatty acid, and phospholipid hydroperoxides. It may play a role in the regulation of phospholipid turnover as well as in protection against oxidative injury (provided by RefSeq, Jul 2008).

Rabbit anti Human peroxiredoxin 6 antibody detects a band of 25 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Rabbit anti peroxiredoxin 6 detects a band of approximately 25 kDa in HeLa cell lysates

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Guarantee 12 months from date of despatch

Acknowledgements PrecisionAb is a trademark of Bio-Rad Laboratories

Health And Safety Information Material Safety Datasheet documentation #20359 available at: Antibody (20359): <https://www.bio-rad-antibodies.com/uploads/MSDS/20359.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M370992:200529'

Printed on 22 Mar 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)