

## Datasheet: VPA00739

<b>Description:</b>	RABBIT ANTI THIOREDOXIN
<b>Specificity:</b>	THIOREDOXIN
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Mouse <b>N.B.</b> Antibody reactivity and working conditions may vary between species.
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Rabbit polyclonal antibody purified by affinity chromatography on Protein A
<b>Buffer Solution</b>	HEPES buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ) 1% Bovine Serum Albumin <50% Glycerol
<b>Immunogen</b>	<i>E. coli</i> -derived recombinant protein of human thioredoxin
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P10599</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">7295</a> TXN <a href="#">Related reagents</a>

<b>Synonyms</b>	TRDX, TRX, TRX1
<b>Specificity</b>	<p><b>Rabbit anti Human thioredoxin antibody</b> recognizes thioredoxin also known as ADF, ATL-derived factor, SASP delta 3, surface-associated sulphhydryl protein or TXN delta 3.</p> <p>The protein encoded by this gene acts as a homodimer and is involved in many redox reactions. The encoded protein is active in the reversible S-nitrosylation of cysteines in certain proteins, which is part of the response to intracellular nitric oxide. This protein is found in the cytoplasm. Two transcript variants encoding different isoforms have been found for this gene (provided by RefSeq, Oct 2011).</p> <p>Rabbit anti Human thioredoxin antibody detects a band of 12 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
<b>Western Blotting</b>	Anti thioredoxin antibody recognizes a band of approximately 12 kDa in HepG2 cell lysates
<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles
<b>Guarantee</b>	12 months from date of despatch
<b>Acknowledgements</b>	PrecisionAb is a trademark of Bio-Rad Laboratories
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10088 available at: Antibody (10088): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10088.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10088.pdf</a>
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

'M371012:200529'

Printed on 11 Aug 2020