

Datasheet: VPA00742

BATCH NUMBER 170504

Description:	RABBIT ANTI SUGT1
Specificity:	SUGT1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by affinity chromatography on immunogen
Buffer Solution	TRIS buffered glycine
Preservative Stabilisers	0.01% Thiomersal <10% Glycerol
Immunogen	Recombinant protein encompassing part of the center region of human SUGT1
External Database Links	UniProt: Q9Y2Z0 Related reagents

Entrez Gene:[10910](#) [SUGT1](#) [Related reagents](#)

Specificity	Rabbit anti Human SUGT1 antibody recognizes suppressor of G2 allele of SKP1 homolog, also known as SGT1B protein or putative 40-6-3 protein. SUGT1 is homologous to the yeast gene SGT1, which encodes a protein involved in kinetochore function and required for the G1/S and G2/M transitions. Two transcript variants encoding different isoforms have been found for SUGT1 (provided by RefSeq, Aug 2008). Rabbit anti Human SUGT1 homolog antibody detects a bands of 38 and 41 kDa and has been extensively validated for western blotting using whole cell lysates.
Western Blotting	Anti SUGT1 antibody recognizes bands of of approximately 38 and 41 kDa in Jurkat cell lysates
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles
Guarantee	12 months from date of despatch
Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories
Health And Safety Information	Material Safety Datasheet documentation #20357 available at: https://www.bio-rad-antibodies.com/SDS/VPA00742 Antibody (20357)
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)**North & South** Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M371015:200529'

Printed on 13 Aug 2023