

Datasheet: VPA00753

Description:	RABBIT ANTI ELF3
Specificity:	ELF3
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by affinity chromatography on immunogen.
Buffer Solution	Phosphate buffered saline.
Preservative	0.01% Thiomersal
Stabilisers	1% Bovine Serum Albumin <20% Glycerol
Immunogen	Synthetic peptide encompassing part of the C terminal region of human ELF3
External Database Links	UniProt: P78545 Related reagents Entrez Gene:

Synonyms	ERT, ESX, JEN
-----------------	---------------

Specificity	Rabbit anti Human ELF3 antibody recognizes the ETS-related transcription factor Elf-3, also known as E74-like factor 3 (ets domain transcription factor), epithelial-restricted with serine box, epithelium-restricted Ets protein ESX or epithelium-specific Ets transcription factor 1. Rabbit anti Human ELF3 antibody detects a band of 47 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.
--------------------	--

Western Blotting	Anti ELF3 antibody recognizes a band of approximately 47 kDa in A431 cell lysates.
-------------------------	--

Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	---

Guarantee	12 months from date of despatch.
------------------	----------------------------------

Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories.
-------------------------	---

Health And Safety Information	Material Safety Datasheet documentation #10096 available at: https://www.bio-rad-antibodies.com/SDS/VPA00753 Antibody (10096)
--------------------------------------	--

Regulatory	For research purposes only.
-------------------	-----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M403367:220721'

Printed on 13 Aug 2023