

## Datasheet: VPA00789KT

<b>Description:</b>	TSPO ANTIBODY WITH CONTROL LYSATE
<b>Specificity:</b>	TSPO
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb™ Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	2 Westerns

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**PrecisionAb antibodies have been extensively [validated for the western blot application](#).** The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Mouse <b>N.B.</b> Antibody reactivity and working conditions may vary between species.
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	20µl Rabbit polyclonal antibody purified by affinity chromatography on immunogen
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ) 2% Sucrose
<b>Immunogen</b>	Synthetic peptide directed towards the C terminal region of human TSPO
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P30536</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">706</a> TSPO <a href="#">Related reagents</a>
<b>Synonyms</b>	BZRP, MBR

**Specificity** **Rabbit anti Human TSPO antibody** recognizes TSPO, also known as translocator protein, benzodiazepine peripheral binding site, mitochondrial benzodiazepine receptor or peripheral-type benzodiazepine receptor.

Present mainly in the mitochondrial compartment of peripheral tissues, the protein encoded by TSPO interacts with some benzodiazepines and has different affinities than its endogenous counterpart. The protein is a key factor in the flow of cholesterol into mitochondria to permit the initiation of steroid hormone synthesis. Alternatively spliced transcript variants have been reported; one of the variants lacks an internal exon and is considered non-coding, and the other variants encode the same protein (provided by RefSeq, Feb 2012).

Rabbit anti Human TSPO antibody detects a band of 18 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting** Anti TSPO antibody recognizes a band of approximately 18 kDa in C2C12 cell lysates

---

**Instructions For Use** Please refer to the [PrecisionAb western blotting protocol](#). For additional information on secondary antibody dilution and exposure time see product web page.

---

**Lysate Composition** 400µg C2C12 lysate lyophilized in RIPA buffer

---

**Lysate Reconstitution** - If using DDT reconstitute the lyophilized lysate with 190µl DI H<sub>2</sub>O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.  
- If using BME reconstitute the lyophilized lysate with 180µl DI H<sub>2</sub>O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.  
Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.

---

**Storage** Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.

---

**Guarantee** As supplied, 12 months from date of despatch

---

**Acknowledgements** PrecisionAb™ is a trademark of Bio-Rad Laboratories.

---

**Health And Safety Information** Material Safety Datasheet documentation #10045 #10561 available at:  
Antibody (10045): <https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf>  
Lysate Material (10561): <https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf>

---

**Regulatory** For research purposes only

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M348009:190214'

