

Datasheet: VPA00825

Description:	RABBIT ANTI PSMD7
Specificity:	PSMD7
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Mouse

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Rabbit polyclonal antibody purified by affinity chromatography on immunogen

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide
<50% Glycerol

Immunogen

Recombinant protein of human PSMD7

External Database**Links****UniProt:**

[P51665](#) [Related reagents](#)

Entrez Gene:

[5713](#) PSMD7 [Related reagents](#)

Synonyms

MOV34L

Specificity

Rabbit anti Human PSMD7 antibody recognizes PSMD7, also known as 26S proteasome non-ATPase regulatory subunit 7, 26S proteasome regulatory subunit S12, 26S proteasome regulatory subunit rpn8, MOV34 and proteasome subunit p40.

The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. PSMD7 encodes a non-ATPase subunit of the 19S regulator. A pseudogene has been identified on chromosome 17 (provided by RefSeq, Jul 2008).

Rabbit anti Human PSMD7 antibody detects a band of 37 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti PSMD7 antibody recognizes a band of approximately 37 kDa in NIH/3T3 cell lysates

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles

Guarantee

12 months from date of despatch

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories

Health And Safety Information

Material Safety Datasheet documentation #10049 available at:
10049: <https://www.bio-rad-antibodies.com/uploads/MSDS/10049.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M380328:210507'

Printed on 21 Mar 2022

