

## Datasheet: VPA00889

<b>Description:</b>	RABBIT ANTI ALPHA ACTININ-1
<b>Specificity:</b>	ALPHA ACTININ-1
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb Polyclonal
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

**The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range.** Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

### Target Species

Human

### Species Cross Reactivity

Reacts with: Mouse

**N.B.** Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

### Product Form

Purified IgG - liquid

### Preparation

Rabbit polyclonal antibody purified by affinity chromatography on immunogen

### Buffer Solution

Phosphate buffered saline

### Preservative Stabilisers

0.09% Sodium Azide  
<50% Glycerol

### Immunogen

Recombinant protein of human alpha actinin-1

---

**External Database****Links****UniProt:**

[P12814](#)   [Related reagents](#)

**Entrez Gene:**

[87](#)   ACTN1   [Related reagents](#)

---

**Specificity**

**Rabbit anti Human alpha actinin-1 antibody** recognizes alpha-actinin-1, also known as F-actin cross-linking protein, actinin 1 smooth muscle and alpha-actinin cytoskeletal isoform.

Alpha actinins belong to the spectrin gene superfamily which represents a diverse group of cytoskeletal proteins, including the alpha and beta spectrins and dystrophins. Alpha actinin is an actin-binding protein with multiple roles in different cell types. In nonmuscle cells, the cytoskeletal isoform is found along microfilament bundles and adherens-type junctions, where it is involved in binding actin to the membrane. In contrast, skeletal, cardiac, and smooth muscle isoforms are localized to the Z-disc and analogous dense bodies, where they help anchor the myofibrillar actin filaments. ACTN1 encodes a nonmuscle, cytoskeletal, alpha actinin isoform and maps to the same site as the structurally similar erythroid beta spectrin gene. Three transcript variants encoding different isoforms have been found for ACTN1 (provided by RefSeq, Jul 2008).

Rabbit anti Human alpha actinin-1 antibody detects a band of 100 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

---

**Western Blotting**

Anti alpha actinin-1 antibody recognizes a band of approximately 100 kDa in HeLa cell lysates

---

**Storage**

Store undiluted at -20°C, avoiding repeated freeze thaw cycles

---

**Guarantee**

12 months from date of despatch

---

**Acknowledgements**

PrecisionAb is a trademark of Bio-Rad Laboratories

---

**Health And Safety Information**

Material Safety Datasheet documentation #20359 available at: Antibody (20359): <https://www.bio-rad-antibodies.com/uploads/MSDS/20359.pdf>

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

**North & South**   Tel: +1 800 265 7376

**America**   Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

**Printed on 22 Mar 2021**

---

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)