

Datasheet: VPA00896

BATCH NUMBER 180704

Description:	RABBIT ANTI OPA1
Specificity:	OPA1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Rabbit polyclonal antibody purified by affinity chromatography on immunogen

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide
<50% Glycerol

Immunogen Recombinant protein of human OPA1

External Database

Links

UniProt:

[O60313](#)

[Related reagents](#)

Entrez Gene:

[4976](#)

OPA1

[Related reagents](#)

Synonyms

KIAA0567

Specificity

Rabbit anti Human OPA1 antibody recognizes OPA1, also known as dynamin-like 120 kDa protein-mitochondrial, dynamin-like guanosine triphosphatase, mitochondrial dynamin-like GTPase, NTG and optic atrophy protein 1.

OPA1 product is a nuclear-encoded mitochondrial protein with similarity to dynamin-related GTPases. It is a component of the mitochondrial network. Mutations in OPA1 have been associated with optic atrophy type 1, which is a dominantly inherited optic neuropathy resulting in progressive loss of visual acuity, leading in many cases to legal blindness. Multiple transcript variants encoding different isoforms have been found for OPA1 (provided by RefSeq, Mar 2009).

Rabbit anti Human OPA1 antibody detects bands of 80-100 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting

Anti OPA1 antibody recognizes a band of approximately 80-100 kDa in LNCaP cell lysates

Storage

Store undiluted at -20°C, avoiding repeated freeze thaw cycles

Guarantee

12 months from date of despatch

Acknowledgements

PrecisionAb is a trademark of Bio-Rad Laboratories

Health And Safety Information

Material Safety Datasheet documentation #10049 available at:
<https://www.bio-rad-antibodies.com/SDS/VPA00896>
10049

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M380297:210507'

