

Datasheet: VPA00943

BATCH NUMBER 180704

Description:	RABBIT ANTI PSIP1
Specificity:	PSIP1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by affinity chromatography on immunogen
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Stabilisers	<50% Glycerol
Immunogen	Recombinant protein of human PSIP1

External Database Links

UniProt:

[O75475](#)

[Related reagents](#)

Entrez Gene:

[11168](#) PSIP1 [Related reagents](#)

Synonyms	DFS70, LEDGF, PSIP2
-----------------	---------------------

Specificity	Rabbit anti Human PSIP1 antibody recognizes PSIP1, also known as PC4 & SFRS1 interacting protein 1, LEDGF and DFS70. Rabbit anti Human PSIP1 detects a band of 78 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.
--------------------	---

Western Blotting	Anti PSIP1 antibody recognizes a band of approximately 78 kDa in A172 cell lysates
-------------------------	--

Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles
----------------	--

Guarantee	12 months from date of despatch
------------------	---------------------------------

Acknowledgements	PrecisionAb is a trademark of Bio-Rad Laboratories
-------------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10049 available at: https://www.bio-rad-antibodies.com/SDS/VPA00943 10049
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M380310:210507'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)