

Datasheet: VPA00945

Description:	RABBIT ANTI SND1
Specificity:	SND1
Format:	Purified
Product Type:	PrecisionAb Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

The PrecisionAb label is reserved for antibodies that meet the defined performance criteria within Bio-Rad's ongoing antibody validation programme. Click [here](#) to learn how we validate our PrecisionAb range. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependent on sample type.

Target Species

Human

Species Cross Reactivity

Reacts with: Rat

N.B. Antibody reactivity and working conditions may vary between species. Cross reactivity is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information.

Product Form

Purified IgG - liquid

Preparation

Rabbit polyclonal antibody purified by affinity chromatography on immunogen

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.09% Sodium Azide
<50% Glycerol

Approx. Protein

IgG concentration 0.5 mg/ml

Concentrations

Immunogen Recombinant protein of human SND1

External Database Links

UniProt:

[Q7KZF4](#) [Related reagents](#)

Entrez Gene:

[27044](#) SND1 [Related reagents](#)

Synonyms TDRD11

Specificity **Rabbit anti Human SND1 antibody** recognizes the Staphylococcal nuclease domain-containing protein 1, also known as EBNA2 coactivator p100, SND1-BRAF fusion, Tudor-SN and tudor domain-containing protein 11.

Rabbit anti Human SND1 antibody detects a band of 120 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.

Western Blotting Anti SND1 antibody recognizes a band of approximately 120 kDa in HepG2 cell lysates

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles

Guarantee 12 months from date of despatch

Acknowledgements PrecisionAb is a trademark of Bio-Rad Laboratories

Health And Safety Information Material Safety Datasheet documentation #10049 available at: <https://www.bio-rad-antibodies.com/SDS/VPA00945>
10049

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M429578:240410'

Printed on 10 Apr 2024
